

# Arterial Spin Labeling

Bradley MacIntosh and Guocheng Jiang

# Learning Objectives


- ASL has been around for over two-decades and is still going strong
- Non-invasive and quantitative, but unique ASL-CBF contrast
  - ASL using a 3 Tesla MRI scanner is the most common application
- There are several ASL image analysis software options, but we are biased to using tools from the FMRIB Software Library
- Showcase diverse applications of the ASL researches across the lifespan and diseases.

# Lecture 2 – How does ASL perfusion MRI work?

Bradley J MacIntosh, PhD  
Senior Scientist, Sunnybrook Research Institute  
Professor, Medical Biophysics, University of Toronto

# ASL relies on a few principles to give perfusion contrast


Spin tagging



Out of volume preparation


<https://mriquestions.com/motsa.html>

Time of flight MR angiography



# ASL relies on a few principles to give perfusion contrast

Rapid multi-slice  
imaging




Functional MRI

Diffusion MRI


Echo planar  
imaging

# ASL relies on a few principles to give perfusion contrast

Spin tagging



Rapid multi-slice imaging



Control and Tag pairs




Out of volume preparation

<https://mriquestions.com/motsa.html>

Echo planar imaging

Serial imaging

## Three stages of a typical ASL experiment: Labeling, Waiting, and Imaging






★ In **continuous ASL (cASL)** and **pseudo-continuous ASL (pcASL)**, a continuous RF pulse will be applied, and proton magnetization will be inverted as blood travel through the labeling plane

In **pulsed ASL (pASL)**, a short RF pulse (~10 ms) is used to invert all the magnetization of blood water in the labeling plane.



## ASL is a non-invasive measurement of cerebral blood perfusion



- **Spin:** Hydrogen (proton) nuclei.
- **Arterial spin labeling:** To (magnetically) label the hydrogen nuclei in the arterial blood.
- An RF field will alter the magnetization of hydrogen nuclei within water that pass through the labeling plane (neck).





# Theory, principles, and applications of ASL date back 2 decades

## *Two of the classic ASL papers*

1996 JCBFM

### Reduced Transit-Time Sensitivity in Noninvasive Magnetic Resonance Imaging of Human Cerebral Blood Flow

\*D. C. Alsop and \*†J. A. Detre



**FIG. 1.** Schematic diagrams of the continuous arterial tagging perfusion imaging models. **A:** In previous implementations, arterial spins are inverted as they pass the inversion plane. Tagged spins gradually relax back towards equilibrium as they diffuse into the tissue.  $\delta$  is defined as the transit time from the inversion plane to the tissue. **B:** When a delay of duration,  $w$ , is introduced, the proximal artery signal is eliminated, and uninverted spins may enter the tissues before imaging. If  $\delta$  increases, fewer uninverted spins arrive in the tissue resulting in an increase in signal. This increase largely cancels the loss of signal resulting from increased  $T_1$  decay of tagged spins before they arrive in the tissue. **C:** In a more complex model, the contribution of intraluminal spins is explicitly accounted for by separating the "tissue" seen in (a) and (b) into a vascular compartment and a true tissue compartment. A second transit time,  $\delta_a$ , from the inversion plane to the arterial vascular compartment, is defined.


# Theory, principles, and applications of ASL date back 2 decades

## *Two of the classic ASL papers*


1998 Magnetic Resonance in Medicine

### **A General Kinetic Model for Quantitative Perfusion Imaging with Arterial Spin Labeling**

Richard B. Buxton, Lawrence R. Frank, Eric C. Wong, Bettina Siewert, Steven Warach, Robert R. Edelman



### ASL “kinetic model”




# ASL remains primarily a research tool

2006 Stroke

## Comparative Overview of Brain Perfusion Imaging Techniques

Max Wintermark, MD; Musa Sesay, MD; Emmanuel Barbier, PhD; Katalin Borbély, MD, PhD;  
William P. Dillon, MD; James D. Eastwood, MD; Thomas C. Glenn, MD; Cécile B. Grandin, MD, PhD;  
Salvador Pedraza, MD; Jean-François Soustiel, MD; Tadashi Nariai, MD, PhD; Greg Zaharchuk, MD, PhD;  
Jean-Marie Caillé, MD; Vincent Dousset, MD; Howard Yonas, MD



PET is  
the gold  
standard

TABLE 1. Overview of the Imaging Techniques Dedicated to Brain Hemodynamics

	Brain Perfusion Imaging Techniques						
	PET	SPECT	XeCT	PCT	DSC	ASL	Doppler
Feasibility	Adults (and children for static exams)	Ionizing radiation				Injection	
Age range	Adults (and children for static exams)	Adults (and children)	Adults (and children)	Adults (and children)	Adults (and children)	Adults + children	Adults + children
Bedside	No	In some instances	No	No	No	No	Yes
Contrast material	$^{15}\text{O}_2$ , $\text{C}^{15}\text{O}_2$ , $\text{H}_2^{15}\text{O}$	$^{133}\text{Xe}$ , $^{99\text{m}}\text{Tc-HMPAO}$ , $^{99\text{m}}\text{Tc-ECD}$ , $^{123}\text{I-IMP}$ (diffusible)	Stable xenon gas (diffusible)	Iodinated contrast material (nondiffusible)	gadolinium chelate (nondiffusible)	None (endogenous contrast)	None (endogenous contrast)
Radiation/study	0.5–2 mSv	3.5–12 mSv	3.5–10 mSv	2–3 mSv	None	None	None
Data acquisition	5–9 min	10–15 min	10 min	40 sec	1 min	5–10 min	10–20 min
Data processing	5–10 min	5 min	10 min	5 min	5 min	5 min	None
Interpretation							
Mathematical model	Kety–Schmidt model	Principle of chemical microspheres for $^{99\text{m}}\text{Tc}$ tracers, Kety–Schmidt model for $^{133}\text{Xe}$ and $^{123}\text{I-IMP}$	Kety–Schmidt model	Meier–Zierler model	Meier–Zierler model	Meier–Zierler model	Other
Assessed parameters	CBF, CBV, rOEF, glucose metabolism	CBF	CBF	CBF, CBV, MTT, TTP, permeability map	CBF, CBV, MTT, TTP, permeability map	CBF	ICA BFV
Large vessels*	No influence on results	No influence on results	No influence on results	Influence results	Influence results	No influence on results	Not applicable
Quantitative accuracy	Yes	Yes for $^{133}\text{Xe}$ and $^{123}\text{I-IMP}$ ; no for the others tracers	Yes	Yes	Not in daily practice	Yes	Yes for hemispheric CBF
Including for low perfused areas†	Yes	Not applicable	Yes	Yes	Not applicable	Not <10 mL/min/100 g	Not applicable
Reproducibility	5%	10%	12%	10–15%	10–15%	10%	5%
Brain coverage	Whole brain	Whole brain	6-cm thickness	4–5 cm thickness	Whole brain	Whole brain	One measurement for each hemisphere
Spatial resolution	4–6 mm	4–6 mm	4 mm	1–2 mm	2 mm	2 mm	Not applicable
Minimal time interval between 2 successive exams	10 min	10 min (split-dose technique for $^{99\text{m}}\text{Tc-HMPAO}$ , $^{99\text{m}}\text{Tc-ECD}$ and $^{123}\text{I-IMP}$ )	20 min	10 min	25 min	0 min	0 min
Clinical applications							
Clinical fields	Chronic cerebrovascular disorders	(Acute and) chronic cerebrovascular disorders	Acute and chronic cerebrovascular disorders	Acute and chronic cerebrovascular disorders	Acute and chronic cerebrovascular disorders	Chronic cerebrovascular disorders	Acute cerebrovascular disorders
		Trauma	Trauma	Trauma	Trauma	Trauma	Trauma
	Dementia and psychiatric diseases	Dementia and psychiatric diseases	Vasospasm	Vasospasm	Vasospasm	Neurodegenerative disorders	Vasospasm
	Epilepsy	Epilepsy	Epilepsy				
	Brain tumors			Brain tumors	Brain tumors	Brain tumors	
	Brain activation studies	Brain activation studies				Brain activation studies	
Emergency setting	No	In some instances	Yes	Yes	Yes	Yes	Yes

# ASL remains primarily a research tool

2006 Stroke

## Comparative Overview of Brain Perfusion Imaging Techniques

Max Wintermark, MD; Musa Sesay, MD; Emmanuel Barbier, PhD; Katalin Borbély, MD, PhD;  
William P. Dillon, MD; James D. Eastwood, MD; Thomas C. Glenn, MD; Cécile B. Grandin, MD, PhD;  
Salvador Pedraza, MD; Jean-François Soustiel, MD; Tadashi Narai, MD, PhD; Greg Zaharchuk, MD, PhD;  
Jean-Marie Caillé, MD; Vincent Dousset, MD; Howard Yonas, MD





TABLE 1. Overview of the Imaging Techniques Dedicated to Brain Hemodynamics

	Brain Perfusion Imaging Techniques				MRI		Doppler
	PET	SPECT	XeCT	PCT	DSC	ASL	
Feasibility					Injection		
Age range	Adults (and children for static exams)	Adults (and children)	Adults (and children)	Adults (and children)	Adults (and children)	Adults+children	Adults+children
Bedside	No	In some instances	No	No	No	No	Yes
Contrast material	$^{15}\text{O}_2$ , $\text{C}^{15}\text{O}_2$ , $\text{H}_2^{15}\text{O}$	$^{133}\text{Xe}$ , $^{99\text{m}}\text{Tc-HMPAO}$ , $^{99\text{m}}\text{Tc-ECD}$ , $^{123}\text{I-IMP}$ (diffusible)	Stable xenon gas (diffusible)	Iodinated contrast material (nondiffusible)	gadolinium chelate (nondiffusible)	None (endogenous contrast)	None (endogenous contrast)
Radiation/study	0.5–2 mSv	3.5–12 mSv	3.5–10 mSv	2–3 mSv	None	None	None
Data acquisition	5–9 min	10–15 min	10 min	40 sec	1 min	5–10 min	10–20 min
Data processing	5–10 min	5 min	10 min	5 min	5 min	5 min	None
Interpretation							
Mathematical model	Kety–Schmidt model	Principle of chemical microspheres for $^{99\text{m}}\text{Tc}$ tracers, Kety–Schmidt model for $^{133}\text{Xe}$ and $^{123}\text{I-IMP}$	Kety–Schmidt model	Meier–Zierler model	Meier–Zierler model	Meier–Zierler model	Other
Assessed parameters	CBV, CBF, rOEF, glucose metabolism	CBF	CBF	CBF, CBV, MTT, TTP, permeability map	CBF, CBV, MTT, TTP, permeability map	CBF	ICA BFV
Large vessels*	No influence on results	No influence on results	No influence on results	Influence results	Influence results	No influence on results	Not applicable
Quantitative accuracy	Yes	Yes for $^{133}\text{Xe}$ and $^{123}\text{I-IMP}$ ; no for the others tracers	Yes	Yes	Not in daily practice	Yes	Yes for hemispheric CBF
Including for low perfused areas†	Yes	Not applicable	Yes	Yes	Not applicable	Not <10 mL/min/100 g	Not applicable
Reproducibility	5%	10%	12%	10–15%	10–15%	10%	5%
Brain coverage	Whole brain	Whole brain	6-cm thickness	4–5 cm thickness	Whole brain	Whole brain	One measurement for each hemisphere
Spatial resolution	4–6 mm	4–6 mm	4 mm	1–2 mm	2 mm	2 mm	Not applicable
Minimal time interval between 2 successive exams	10 min	10 min (split-dose technique for $^{99\text{m}}\text{Tc-HMPAO}$ , $^{99\text{m}}\text{Tc-ECD}$ and $^{123}\text{I-IMP}$ )	20 min	10 min	25 min	0 min	0 min
Clinical applications							
Clinical fields	Chronic cerebrovascular disorders	(Acute and) chronic cerebrovascular disorders	Acute and chronic cerebrovascular disorders	Acute and chronic cerebrovascular disorders	Acute and chronic cerebrovascular disorders	Chronic cerebrovascular disorders	Acute cerebrovascular disorders
		Trauma	Trauma	Trauma	Trauma	Trauma	Trauma
	Dementia and psychiatric diseases	Dementia and psychiatric diseases	Vasospasm	Vasospasm	Vasospasm	Neurodegenerative disorders	Vasospasm
	Epilepsy	Epilepsy	Epilepsy				
	Brain tumors			Brain tumors	Brain tumors	Brain tumors	
	Brain activation studies	Brain activation studies				Brain activation studies	
Emergency setting	No	In some instances	Yes	Yes	Yes	Yes	Yes

# DSC versus ASL: Which MRI perfusion approach is more popular?

3,023 results


RESULTS BY YEAR



<https://pubmed.ncbi.nlm.nih.gov/?term=%22dynamic+susceptibility+contrast%22+or+%28%22DSC%22+and+MRI%29&sort=pubdate>

5,207 results


RESULTS BY YEAR



<https://pubmed.ncbi.nlm.nih.gov/?term=%22arterial+spin+labeling%22+or+%28%22ASL%22+and+MRI%29+or+%22arterial+spin+tagging%22&sort=pubdate>



## Pseudo-continuous ASL (pcASL) imaging was recommended in 2015 and this helped streamline neuroimaging research



### Recommended Implementation of Arterial Spin Labeled Perfusion MRI for Clinical Applications: A consensus of the ISMRM Perfusion Study Group and the European Consortium for ASL in Dementia

David C. Alsop<sup>1,\*</sup>, John A. Detre<sup>2</sup>, Xavier Golay<sup>3</sup>, Matthias Günther<sup>4,5,6</sup>, Jeroen Hendrikse<sup>7</sup>, Luis Hernandez-Garcia<sup>8</sup>, Hanzhang Lu<sup>9</sup>, Bradley J. MacIntosh<sup>10,11</sup>, Laura M. Parkes<sup>12</sup>, Marion Smits<sup>13</sup>, Matthias J. P. van Osch<sup>14</sup>, Danny JJ Wang<sup>15</sup>, Eric C. Wong<sup>16,†</sup>, and Greg Zaharchuk<sup>17</sup>

- **pcASL** can be applied on clinical scanners without special hardware.
- There is a well-defined labeling duration (LD) for perfusion quantification and to maximize the SNR.
- Relatively high labeling efficiency.

ISMRM

Perfusion

STUDY GROUP




# ASL topics not discussed

**Spatial resolution** – can we push towards smaller voxel dimensions?

**Main magnetic field** – what is the status of ASL at 7 Tesla

**Pulse sequences** – time resolved angiography, vessel selective, velocity selective



## So What Is the **ASL-CBF** Signal Anyway?

- ASL-CBF signal depends on water, which is a freely diffusible tracer.
- ASL-CBF signal is a measure of tissue perfusion.
- ASL has its wrinkles that could include arterial transit time and partial volume error effects
  - Visually-inspect your ASL-CBF maps

