

Structural MRI Module

Introduction to MRI physics & hardware

Jason Rock



Sunnybrook
RESEARCH INSTITUTE



Medical Biophysics
UNIVERSITY OF TORONTO


Outline

1. Basic principles of MRI
2. MRI hardware – what's inside the scanner
3. Safety considerations
4. Making an image
5. What's the deal with 7 Tesla?
6. Why we use MRI

Basic Principles of MRI

Your body is full of water

- Specifically, hydrogen atoms
- Hydrogen atoms have a **magnetic moment**




MRI In Practice, Westbrook 5th edition




WATER IN THE HUMAN BODY

Brain	75% Water
Blood	83% Water
Heart	79% Water
Bones	22% Water
Muscles	75% Water
Liver	85% Water
Kidneys	83% Water

Basic Principles of MRI




Random alignment
No external field




Alignment
External magnetic field

MRI In Practice, Westbrook 5th edition


Similar to a compass needle aligning with Earth's magnetic field




Basic Principles of MRI



Random alignment
No external field



Alignment
External magnetic field





MRI In Practice, Westbrook 5th edition

Basic Principles of MRI - Resonance

A resonant radio-frequency (RF) pulse can topple the net magnetization vector (NMV) over

- The NMV will precess in the B_0 field, inducing a signal in the same coil



MRI In Practice, Westbrook 5th edition



MRI Hardware

The MRI has 3 main components

1. Static B_0 field
2. Gradient magnetic fields
3. RF transmission/receive



Siemens




MRI In Practice, Westbrook 5th edition

MRI Hardware – B_0 field

A superconducting coil of wire creates the B_0 field

- Clinical field strengths range from 0.5T to 3T
- The magnet is always on
- Additional conducting coils help shield and shim the magnet





MRI In Practice, Westbrook 5th edition

MRI Hardware – Gradients

The gradient coils encode spatial information in the MR signal

- Variation in field strength in all three axis
- Culprit for the loud noise in MRI




MRI In Practice, Westbrook 5th edition

MRI Hardware – RF transmission/receive

RF antennas are used to transmit the NMV, but also used to receive the induced MR signal

- The closer the coil to the subject, the higher the SNR



Safety Considerations

BBC

Man wearing heavy metallic necklace dies after being sucked into MRI machine

21 July 2025

Madeline Halpert
BBC News, New York

Share  Save 



Forbes

INNOVATION > HEALTHCARE

Lawyer Dies After Shot By His Own Concealed Gun Triggered By MRI Scanner

By Bruce Y. Lee, Senior Contributor. © Bruce Y. Lee, M.D., MBA, covers health... 

[Follow Author](#)

Published Feb 12, 2023 at 09:02am EST, Updated Feb 12, 2023 at 02:07pm EST



Safety Considerations

The MRI is a strong magnet

- We always say “absolute no metals inside an MRI”
- However, some metals are unaffected by magnetic fields

Ferromagnetic Metal

A metal that strongly attracts a magnetic field

- Iron
- Nickel
- Cobalt
- Steel

Non-ferromagnetic Metal

A metal that does not exhibit strong magnetic properties

- Copper
- Aluminum
- Gold
- Titanium




Safety Considerations

However, even non-ferromagnetic metals can be dangerous

- The RF transmission can couple with metals causing immense heating
- Wires and metallic leads are most susceptible to heating

Fork in microwave




Henderson, 2005. Permanent neurological deficit related to magnetic resonance imaging in a patient with implanted deep brain stimulation electrodes for Parkinson's disease: Case report

Safety Considerations

However, even *normal* tissue can heat up

- The Specific Absorption Rate (SAR) is a measure of RF energy absorption
- Calculated by scanner, given sequence parameters, patient body mass & height





<https://itis.swiss/customized-research/mri-safety-evaluation/>

Making an Image

The gradients add spatial localization to the MR signal

- Each tissue has a unique MR signal




MRI In Practice, Westbrook 5th edition

Making an Image


The MR signal will decay due to

- Longitudinal magnetization recovery (**T1**)
- Loss of transverse magnetization coherence (**T2**)

T2 relaxation



T1 relaxation



MRI In Practice, Westbrook 5th edition

Making an Image

Each tissue has a unique T1 and T2 relaxation constant

Tumor Detection by Nuclear Magnetic Resonance

Raymond Damadian

Table 1. Spin-lattice (T_1) and spin-spin (T_2) relaxation times (in seconds) of normal tissues.

Rat No.	Weight (g)	Tissue							
		Rectus muscle		Liver		Stomach	Small intestine	Kidney	Brain
		T_1	T_2	T_1	T_2	T_1	T_1	T_1	T_1
1	156	0.493	0.050	0.286	0.050	0.272	0.280	0.444	0.573
2	150	.548	.050	.322	.060	.214	.225	.503	.573
3	495	.541	.050	.241	.050	.260	.316	.423	.596
4	233	.576 (0.600)*	.070	.306 (0.287)*	.048	.247 (0.159)*	.316 (0.280)*	.541 (0.530)*	.620 (0.614)*
5	255	.531		.300		.360	.150	.489	.612
Mean and standard error									
		0.538 ± 0.015	0.055 ± 0.005	0.293 ± 0.010	0.052 ± 0.003	0.270 ± 0.016	0.257 ± 0.030	0.480 ± 0.026	0.595 ± 0.007

* Spin-lattice relaxation time after the specimen stood overnight at room temperature.

Table 2. Spin-lattice (T_1) and spin-spin (T_2) relaxation times (in seconds) in tumors.

Rat No.	Weight (g)	T_1	T_2
<i>Walker sarcoma</i>			
6	156	0.700	0.100
7	150	.750	.100
8	495	.794 (0.794)*	.100
9	233	.688	
10	255	.750	
Mean and S.E.		0.736 ± 0.022	.100
P		$< .01 \uparrow$	
<i>Novikoff hepatoma</i>			
11	155	0.798	0.120
12	160	.852	.120
13	231	.827	.115
Mean and S.E.		0.826 ± 0.013	0.118 ± 0.002
P		$< .01 \uparrow$	
<i>Fibroadenoma (benign)</i>			
14		0.448	
15		.537	
Mean		.492	

(measured in 1971)

Making an Image – Pulse Sequence

An MR image is created through the precise **orchestration** of the gradients and RF transmitter in the pulse sequence


RF Transmitter

X-Gradient

Y-Gradient

Z-Gradient

RF Receiver



Time



What's the deal with 7 Tesla?


Higher magnetic field strength allows for

- Higher SNR, CNR, changes in T1 and T2 relaxation constants


7 Tesla MRI allows us to unlock

- The tiniest structures and vessels in the brain
- Subtle tissue differences


3 Tesla




7 Tesla



globus pallidus sub-segments




What's the deal with 7 Tesla?




Rooney WD et al. Magnetic field and tissue dependencies of human brain longitudinal ¹H₂O relaxation in vivo. *Magn Reson Med*. 2007 Feb;57(2):308-18.


3 Tesla



7 Tesla



globus pallidus sub-segments




Why Choose MRI?

MR imaging has exquisite soft tissue contrast


- Gray and white matter have nearly the same density
- CT image results in poor contrast

No ionizing radiation

- Unlimited number of scanning opportunities




Why Choose MRI?




Conventional Imaging
(1972)





Angiogram
(Dumoulin and Hart 1986)




Spectroscopy (Carr 1952)



Diffusion Image
(Bihan 1986)



fMRI
(Seiji Ogawa 1990)



Phase contrast
(~1980s?)